# Empress® System

## Clinical Guide

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# Table of contents

<b>1. Indications</b>	<b>14</b>
<b>2. Clinical procedure</b>	<b>15</b>
Shade selection	15
Preparation	16
Preparation techniques for individual indications:	17
- Veneers	
- Inlays	
- Onlays	
- Crowns	
Die shade selection	19
Impression taking	19
Temporary restoration	19
Cementation	20

# 1. Indications

Indications			
	IPS Empress	IPS Empress 2	
ANTERIOR REGION	Veneers Crowns Bridges	• • —	— • <sup>1)</sup> • <sup>2)</sup>
POSTERIOR REGION	Inlays/Onlays Crowns Bridges	• • —	— • <sup>1)</sup> • <sup>3)</sup>
<p>1) conventional cementation also feasible 2) conventional cementation also feasible; bridges with one pontic 3) conventional cementation also feasible; bridges with one pontic and up to the 2<sup>nd</sup> premolar as the distal abutment</p>			

## 2. Clinical procedure

Starting situation

### Shade selection

If necessary, conduct professional tooth cleaning prior to selecting the shade. The shade of the unprepared tooth is determined with the help of the Chromascop shade guide. Do not overdry the tooth.

Make sure to observe individual characteristics when selecting the shade. If a crown preparation is planned, the incisal and cervical shades should also be determined. In order to achieve optimum shade match, patients should not wear brightly coloured clothes or lipstick during shade selection. The shade should be determined under natural lighting conditions, whenever possible.

Shade determination on the patient

## Preparation

Remove existing restorations and/or caries.

### Cavity preparation

All-ceramic restorations like IPS Empress® and IPS Empress® 2 are fabricated according to the principles of adhesive cementation. To ensure quality restorations, the following principles need to be observed. These points provide an overview of the preparation requirements for indirect restorations seated using the adhesive technique. For further information, we recommend consulting the literature or attending suitable continuing education courses.

### Basic rules:

- Complete a shoulder/chamfer preparation with rounded angles; the width of the circular shoulder/chamfer should be approx. 1 mm.
- Do not prepare feather edges; Use butt joints.
- Block out undercuts.
- Avoid sharp edges and angles to prevent stress and facilitate insertion.
- Strictly observe the necessary minimum wall thickness in order to ensure suitable stability of the restoration.
- Complete the prep with even reduction and maintain the anatomical contours.

Completed inlay preparation

## Preparation techniques for individual indications

### Veneers

If possible, the preparation should be entirely located in the enamel.



The chamfer preparation should involve a clearly defined margin. Two techniques are available. The first involving a wraparound technique onto the lingual surface and the other stopping at the top of the incisal edge. The incisal margins should not be located in the occlusal contact areas.

The dimension of the incisal reduction depends on the desired degree of translucency of the incisal area to be restored. The more transparent the incisal edge of the veneer should appear, the more reduction is required. 1.0-mm to 1.5-mm incisal reduction is recommended.

By providing depth grooves using the appropriate burs, controlled enamel reduction can be achieved.

The minimum preparation thickness is approx. 0.6–1.0 mm depending on the fabrication technique chosen.

Eliminating proximal contacts is not required. Discoloured teeth may require more extensive reduction to allow for adequate masking. Observe an angle of approx. 10°–30° for cervical chamfers or shoulder preparations (same as for crowns).




### Inlays

Avoid occlusal contacts at the margins.

Observe a preparation depth of at least 1.5 mm in the fissure area. The isthmus should be at least 1.5 mm wide.

Complete the proximal box with a slightly diverging angle. Prepare an angle of >60° between the proximal cavity walls and the prospective proximal surfaces of the inlay. Round out internal angles. Cavosurface margins should be prepared as butt joints.



## Onlays

Same procedure as for inlays.

Provide 2 mm space in the area of the cusps and occlusal contacts.


Prepare the shoulder with an incline (10°–30°), to improve the aesthetic appearance of the transitional area between ceramic material and tooth. Onlays are indicated if the preparation margin runs less than approx. 0.5 mm away from the cusp tip, or if the enamel is severely undermined.

## Crowns


Evenly reduce the anatomical shape observing the stipulated minimum wall thickness.

Prepare a chamfer/shoulder preparation (10°–30°) with rounded inner angles. The width of the circular shoulder/chamfer should be approx. 1.0 mm.


Reduce the axial walls of the tooth on the incisal third of the crown by approx. 1.5 mm.

The incisal and/or occlusal reduction should be approx. 1.5–2.0 mm.

For anterior crowns, the labial or palatal/lingual reduction should be approx. 1.0–1.5 mm.

## Die shade selection

The shade of the prepared die may be determined using the die material shade guide (Die Material Shade Guide). With the appropriate shade, the dental technician is able to fabricate aesthetically appealing, true-to-nature restorations. In order to achieve optimum results, the dentin should be moist when the shade is determined.

## Temporary restoration

Inlay and onlay cavities may be restored with a temporary restorative material (e.g. Systemp®.inlay/Systemp®.onlay).

Composite resin temporaries are the best solution for crown and bridge preparations.


Use a temporary eugenol-free cement (e.g. Provilink®) for provisional cementation. Do not use cements that contain eugenol, since this substance may inhibit the polymerization of the luting cement for the permanent restoration.

Temporary restoration fabricated of Systemp®.inlay

## Impression taking


An addition-curing silicone or other suitable impression material may be used for the impression.

Impression of the completed preparation

## Cementation

Generally, adhesive cementation (e.g. using the Variolink® II System) is recommended for metal-free restorations. If this technique is not possible for clinical reasons, crowns and bridges fabricated of IPS Empress 2 may also be seated using the conventional cementation technique in conjunction with glass ionomer cements (e.g. Vivaglass® CEM).

For further information, please refer to the Variolink II Instructions for Use

Incorporated restoration (IPS Empress<sup>®</sup>)

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
**Date information prepared:** 10/2002

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty and are not binding.

Printed in Austria  
© Ivoclar Vivadent AG, Schaan / Liechtenstein  
561526/1002/d/e/3/VVA

The logo consists of a series of small green dots arranged in a stylized, horizontal pattern that tapers to the right. Below this graphic, the brand name "ivoclar vivadent" is written in a lowercase, sans-serif font. The word "ivoclar" is in blue, and "vivadent" is in black. A thin horizontal line extends from the end of "vivadent" to the right. Underneath "vivadent", the word "technical" is written in a smaller, italicized, lowercase sans-serif font.

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**Ivoclar Vivadent S.A.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent, Inc.**